

LSR II Tutorial Checklist

Working in the flow unit:

- No food or beverages are allowed except for closed bottles
- The workstation should be left clean and organized
- Ensure you take your belongings, such as articles, protocols, experiment layout, etc.
- Instructions on the reservation system and working with the KIOSK
- REGISTRATION FORM, Reservation and registration, Safety Instructions, Shadow protocol.

Kiosk

- ☐ Sign in before starting work (verify budget)
- ☐ Sign out after finishing work
- ☐ If you have a tutorial, sign up for both sessions
- ☐ Verify that the next user is coming (after working hours)


Machine

- ☐ Turn on/off procedure
- ☐ Cleaning machine
- ☐ Read sample (overwrite and append)
- ☐ UV laser and HTS (optional)
- ☐ Turn off LSR if the next user has more than 2 hours until their session

Diva

- ☐ Experiment duplication
- ☐ Area, Width, and Height
- ☐ Adjust voltage
- ☐ Graphs and Gates
- ☐ Gates
- ☐ Record mode (storage gate and stopping gate)
- ☐ Compensation (optional)
- ☐ Export data and access data (Responsibility of the data – by the user !!!)


Figure A-10 Anatomy of a pulse



Pulse Measurements

The pulse processors measure pulses by three characteristics: height, area, and width.

Figure A-11 Pulse measurements



- Pulse height is the maximum digitized intensity measured for the pulse.
- Pulse area is an integration of the digitized measures over time.
- Pulse width calculates: $\frac{\text{area}}{\text{height}} \times 64,000$