

## ImageStream Tutorial Checklist

### Working in the flow unit:

- No food or beverages are allowed except for closed bottles
- The workstation should be left clean and organized
- Ensure you take your belongings, such as articles, protocols, experiment layout, etc.
- Instructions on the reservation system and working with the KIOSK
- REGISTRATION FORM, Reservation and registration, Safety Instructions, Shadow protocol.

### Kiosk

- ☐ Sign in before starting work (verify budget)
- ☐ Sign out after finishing work
- ☐ If you have a tutorial, sign up for both sessions
- ☐ Turn off the **stream** if the next user has more than 2 hours until their session
- ☐ Verify that the next user is coming (after working hours)

### Machine

- ☐ Fluids: Sterilizer (0.5% bleach), Cleanser (coulter clenx), Debbubler (70% isopropyl alcohol)  
Rinse (DDW), Sheath (PBSX1), Beads (SpeedBead)
- ☐ Waste –make sure it's not full before turning on the machine
- ☐ Turn on/off procedure (server, ImageStream, pc, Inspire)
- ☐ Cleaning and calibrating the machine
- ☐ If the staff is available, we can power on the machine and run calibration, it's the researcher  
Responsibility to notify us before coming to the FACS unit!

### Inspire


- ☐ Experiment duplication (template)
- ☐ Sample handling (volume, load, and return), Acquisition (counting and recording)
- ☐ Internal beads exclusion
- ☐ Channel to record and eliminate, Brightfield
- ☐ Compensation– Mandatory - the ImageStream is co-linear

- ☐ Adjust laser voltages
- ☐ Magnification (X20, X40, X60)
- ☐ Graphs -Area, Intensity, Aspect ratio, Gradient RMS
- ☐ Gates (creating, renaming, modifying)
- ☐ Record mode (storage gate and stopping gate)
- ☐ Backup data, deleting files on the desktop
- ☐ Optional for big experiment copy to external HDD (Deepfreeze)
- ☐ Cheat sheet - explanation


#### Analysis – IDEAS (schedule a tutorial or time with the staff)

- ☐ File types and usage
- ☐ Wizards
- ☐ Export data and access data
- ☐ Mask and Features
- ☐ Help file (mainly for analysis)


### Basic gating strategy (cells)



Internal beads exclusion



Doublets discrimination




Then continue gating according to your desired dyes and fluorochrome

### ImageStram cheat sheet


\*For compensation turn OFF brightfield and SSC, a prefix of \_noBF will be added to your files

### Basic gating strategy (Beads)




Area of Brightfield

Internal beads exclusion



Area of Brightfield


Doublets discrimination



Gradient RMS of Brightfield

Then continue gating according to your desired dyes and fluorochrome


**Figure A-10** Anatomy of a pulse



## Pulse Measurements

The pulse processors measure pulses by three characteristics: height, area, and width.

**Figure A-11** Pulse measurements



- Pulse height is the maximum digitized intensity measured for the pulse.
- Pulse area is an integration of the digitized measures over time.
- Pulse width calculates:  $\frac{\text{area}}{\text{height}} \times 64,000$