



Date modified: October 2025

Date of the tutorial: _____

Aurora Tutorial Checklist

Working in the flow unit:

- No food or beverages are allowed except for closed bottles
- The workstation should be left clean and organized
- Ensure you take your belongings, such as articles, protocols, experiment layout, etc.
- Instructions on the reservation system and working with the KIOSK
- REGISTRATION FORM, Reservation and registration, Safety Instructions, Shadow protocol.

Kiosk

- ☐ Sign in before starting work (verify budget)
- ☐ Sign out after finishing work
- ☐ If you have a tutorial, sign up for both sessions
- ☐ Turn off Aurora if the next user has more than 4 hours until their session
- ☐ Verify that the next user is coming (after working hours)

Machine

- ☐ Fluidics: Need to wash at least one hour before starting work
- ☐ Run Spectroflo QC Beads daily


Spectroflo

- ☐ Create a new user
- ☐ Experiment duplication (with/without reference/data, template)
- ☐ Sample Acquisition (counting and recording, moving between samples)
- ☐ Unmixing
- ☐ Adjust detector gain (add height and width)
- ☐ Graphs
- ☐ Gates (creating, renaming, modifying)
- ☐ Record mode (storage gate and stopping gate)
- ☐ Export and access data
- ☐ Spectral signatures



- ☐ Help file (LS&E site), full spectrum viewer tool
- ☐ Analysis on spectroflo (analysis station, no user name and password)
Or Flowjo, FCSEXpress and Cytobank
- ☐ Shutdown


Figure A-10 Anatomy of a pulse



Pulse Measurements

The pulse processors measure pulses by three characteristics: height, area, and width.

Figure A-11 Pulse measurements



- Pulse height is the maximum digitized intensity measured for the pulse.
- Pulse area is an integration of the digitized measures over time.
- Pulse width calculates: $\frac{\text{area}}{\text{height}} \times 64,000$